


## Side panels for medicine cabinets

The ISO modular system plastic side panels allow opening the baskets on a 5° as well as a 30° angle, which significantly improves access to supplies in the upper baskets. Stoppers in basket corners ensure the reliable fixing of baskets in an open position.

When storing heavy accessories, the side panels can be equipped with rollers that assist pulling out the baskets.

Thanks to the unique construction, the plastic side panels can be used in passage cabinets where there must be a possibility to pull out the baskets on both sides.




Art. 7900030 Side panel 300, left, dimensions(WxH) 263x640 mm.  
Material ABS, grey.

Art. 7900031 Side panel 300, right, dimensions(WxH) 263x640 mm.  
Material ABS, grey.



Art. 7900032 Side panel 400, dimensions(WxH) 330x640 mm.  
Material ABS, grey. Can be used in pass through cabinets.



Art. 7900033 Side panel 600, dimensions(WxH) 530x640 mm.  
Material ABS, grey. Can be used in pass through cabinets.



Art. 90.0013R Side panel inclined, right, dimensions(WxH) 365x395 mm.

Art. 90.0013L Side panel inclined, left, dimensions(WxH) 365x395 mm.